



ORDERING INFORMATION

TOLL FREE

800. 557. 3587

FAX

720. 873. 0207

AlloMend® Acellular Dermal Matrix

STERILE	PRODUCT DESCRIPTION		
73583008	2 x 4 cm	(8 cm ²)	0.4 - 1.0 mm (M)
73583016	4 x 4 cm	(16 cm ²)	0.4 - 1.0 mm (M)
73583032	4 x 8 cm	(32 cm ²)	0.4 - 1.0 mm (M)
73583128	8 x 16 cm	(128 cm ²), 1:1 meshed	0.4 - 1.0 mm (M)
73083008	2 x 4 cm	(8 cm ²)	1.0 - 2.0 mm (T)
73083024	2 x 12 cm	(24 cm ²)	1.0 - 2.0 mm (T)
73083016	4 x 4 cm	(16 cm ²)	1.0 - 2.0 mm (T)
73083032	4 x 8 cm	(32 cm ²)	1.0 - 2.0 mm (T)
73083048	4 x 12 cm	(48 cm ²)	1.0 - 2.0 mm (T)
73083064	4 x 16 cm	(64 cm ²)	1.0 - 2.0 mm (T)
73083072	6 x 12 cm	(72 cm ²)	1.0 - 2.0 mm (T)
73083096	6 x 16 cm	(96 cm ²)	1.0 - 2.0 mm (T)
73083128	8 x 16 cm	(128 cm ²)	1.0 - 2.0 mm (T)
73183016	4 x 4 cm	(16 cm ²)	2.0 - 3.3 mm (XT)
73183032	4 x 8 cm	(32 cm ²)	2.0 - 3.3 mm (XT)
73183064	4 x 16 cm	(64 cm ²)	2.0 - 3.3 mm (XT)
73183128	8 x 16 cm	(128 cm ²)	2.0 - 3.3 mm (XT)

(M) Medium, (T) Thick, (XT) Extra-Thick

Other sizes and meshing options may be available on request



6278 S Troy Cir
Centennial, CO 80111

For more information,
please call 720. 873. 0213
or visit allosource.org


AlloSource, a non-profit organization, offers more than 200 types of precise bone, skin, soft-tissue and custom-machined allografts for use in an array of life-saving and life-enhancing medical procedures. Committed to honoring the gift of donation, the company delivers unparalleled expertise and customer service to its network of surgeons, partners and the country's most reputable Organ Procurement Organizations.

Leader in cellular tissue processing including fresh skin allografts for severe burns, fresh cartilage tissue for joint repair and adult mesenchymal stem cells.



ALLOMEND[®]

DOING MORE
FOR SOFT TISSUE
REPAIR AND
RECONSTRUCTION.



ACELLULAR DERMAL MATRIX



→ INTRODUCING ALLOMEND®

AlloSource, a non-profit supplier of allografts and one of the nation's leading tissue banks, introduces AlloMend Acellular Dermal Matrix (ADM) bringing the benefits of regenerative medicine to more patients.

Human acellular matrices are used in a broad range of surgical procedures including:

- Hernia repair¹
- Abdominal wall reconstruction¹
- Pelvic organ prolapse¹
- Breast reconstruction²
- Tendon augmentation³
- Rotator cuff repair⁴

→ ACELLULAR REGENERATION

Through a proprietary process, viable cells and cellular elements that are capable of triggering an immunogenic response are removed from donated human dermal tissue, leaving a collagen elastin matrix behind. Upon transplantation, the body's own cells infiltrate and repopulate this three-dimensional scaffold to begin the revascularization and remodeling processes. Acellular allograft matrices, unlike synthetic materials or xenografts, are recognized as human tissue by the body for graft incorporation by the recipient, minimizing the risk of inflammation⁵ or rejection⁶.




→ ALLOMEND BASICS

THOROUGH DECELLULARIZATION

AlloMend ADM utilizes a dynamic tissue cleansing process, without the use of detergents or enzymes. The result is thorough decellularization, but with no harmful residuals in the tissue.

AlloMend processing results in significant removal of cellular debris, including, DNA, RNA, proteins and antigens, without altering the morphological collagen structure.




H&E (hematoxylin and eosin) stain review of "before and after" decellularization process.

HIGH TENSILE STRENGTH

AlloMend ADM exceeds the tensile strength of leading acellular dermal matrices for more assurance in surgical repair of integumental tissue.

Ultimate tensile strength is a standard testing methodology to measure the force needed to stretch and break a biomaterial.



A CLOSER LOOK AT ALLOMEND

- **Flexible and pliable material** - optimal handling characteristics enable precision placement
- **Available in a variety of thicknesses and sizes** - wide range of surgical applications
- **Terminally sterilized to a Sterility Assurance Level (SAL) of 10^{-6} , with e-beam technology** - minimizes infection risk, while avoiding damaging tissue
- **Two-year shelf life in room-temperature conditions** - no special handling or storage required
- **Dual layer packaging** - inner sterile pouch enables placement in sterile OR field
- **Prehydrated** - immediately ready to use, no need to wait for product to rehydrate

SAFETY AND QUALITY OF ALLOMEND

AlloSource requires a comprehensive donor physical assessment and a complete medical and social history to identify and eliminate donors that may be at risk of transmitting certain viruses and diseases. Our donor acceptance criteria are based on regulations established by the U.S. Food and Drug Administration (FDA), the American Association of Tissue Banks (AATB) Standards, as well as additional requirements set by the AlloSource Medical Advisory Board.

Donors must test negative or non-reactive in the following assays:

- | | |
|--|---|
| - Antibody to Hepatitis C (HCV) | - Hepatitis C Virus (HCV NAT) |
| - Antibody to Human Immunodeficiency Virus 1 & 2 (HIV 1 & 2) | - Human Immunodeficiency Virus Type 1 (HIV-1 NAT) |
| - Hepatitis B Core IgG/IgM Antibody (HBcAb) | - Rapid Plasma Reagin or Serologic Test for Syphilis (RPR or STS) |
| - Hepatitis B Surface Antigen (HBsAg) | |

Allograft tissue supplied exclusively by our partner organ procurement organizations (OPOs), located domestically.

1. Pappas G et al. Biological mesh in hernia repair, abdominal wall defects reconstruction and treatment of pelvic organ prolapse: A review of the clinical evidence. *The American Surgeon*. 76(11): 1290-99 (2010).
2. Kokac E, et al. Biologic matrices in oncologic breast reconstruction after mastectomy. *Expert Review of Medical Devices*. 11(1): 65-75 (2014).
3. Wilkins, R. Acellular dermal grafts augmentation in quadriceps tendon rupture repair. *Current Orthopaedic Practice*. 21(3): 315-19 (2010).
4. Barber FA, et al. A prospective, randomized evaluation of acellular human dermal matrix augmentation for arthroscopic rotator cuff repair. *Arthroscopy*. 28(1): 8-15 (2012).
5. Richters C, et al. Development of a dermal matrix from glycerol preserved allogenic skin. *Cell and Tissue Banking*. 9(4): 309-15 (2008).
6. Michael TE. Xenograft risks: What you and your patients need to know. *American Academy of Orthopaedic Surgeons*. www.aaos.org/news/aaosnow/jun09/research3.asp



Sport Medicine Grafts



Frozen

Reference	Description Tendons	Auxiliary
AS10017000	Achilles with Bone block	
AS17917000	Patellar Lgmt. Whole (BTB)	
AS17817000	Patellar Lgmt. Hemi (BTB)	Ligament length > 4.2 cm
AS18017000	Patellar Lgmt. Preshaped (BTB)	Ligament length > 4.2 cm
1964-14	Patellar Lgmt. Preshaped (BTB) / Dowels 10 mm	Ligament length > 4.2 cm
3398-14	Patellar Lgmt. / LARGE 12 mm / Preshaped (BTB) / Dowels 10 mm	Ligament length 3–5 cm, Width 1.2 cm
AS43917000	Peroneus Longus	L = 28–30 cm
1170-14	Peroneus Longus	L = 28–30 cm
AS43917000XL	Peroneus Longus XLarge	L = > 30 cm
1170-14XL	Peroneus Longus XLarge	L = > 30 cm
AS44317002	Gracilis, Single Strand	L = 18–23 cm
1009-14	Gracilis, Single Strand	L = 18–23 cm
AS44317003	Semitendon, Single Strand	L = 21–23 cm
1824-14	Semitendon, Single Strand	L = 21–23 cm
AS44317000	Anterior Tibialis, Single Str.	L = 18–23 cm
AS44317001	Posterior Tibialis, Single Str.	L = 18–23 cm
1162-14R	Patellar Lig w/Extensor (Right)	
1162-14L	Patellar Lig w/Extensor (Left)	
1903-14	Quadriceps w/Bone Block 2.5 cm	L = 21–24 cm

Patellar Ligament - Hemi



Patellar Ligament - Whole



Patellar Ligament - Pre-shaped



Achilles Tendon with Bone Block



Posterior Tibialis



Peroneus Longus



Semi-T



Gracilis



Patellar Ligament with Extensor



Tendon-Quadriceps with Bone Block





Traditional Grafts

TriCortical Ilium Wedges



TriCortical Ilium Wedges

Freeze-Dried	Frozen	Product Description
30515008	30517008	≤ 8 mm (W)
30615013	30617013	9 - 13 mm (W)
30715016	30717016	14 - 16 mm (W)
30815019	30817019	17 - 19 mm (W)
30915022	30917022	20 - 22 mm (W)
31015023	31017023	≥ 23 mm (W)

BiCortical Ilium Strip



TriCortical Ilium Strips

Freeze-Dried	Frozen	Product Description
20131500	20317000	TriCortical Ilium Strip

BiCortical Block

Freeze-Dried	Frozen	Product Description
10615040		≤ 4 cm (L)

Dowel



BiCortical Ilium Strips

Freeze-Dried	Frozen	Product Description
10815041		≥ 4 cm (L)

Dowels w/o Cartilage - UniCortical

Freeze-Dried	Frozen	Product Description
14215012		12 mm
14215013		13 mm
14215014		14 mm
14215015		15 mm
14215018		18 mm

Patella Wedges



Patella Wedges

Freeze-Dried	Frozen	Product Description
41715016		Wedge ≤ 16 mm
41715017		Wedge ≥ 16 mm

Cortical Strut



Cortical Struts

Freeze-Dried	Frozen	Product Description
12115010	12117010	1 cm x ≤ 10 cm, approximate
12115024	12117024	2 cm x 11 - 24 cm, approximate
12115025	12117025	2 cm x ≥ 25 cm, approximate